



## Binocular Digital Stereo Microscope biological Microscopes 5x / 10x Rohs A22.1209

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1209
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 5-80X
- Applications: Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,  
stereo binocular microscopes**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Binocular Digital Stereo Microscope biological Microscopes Rohs A22.1209

#### Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

#### Specifications:

A22.1209 Specification									
Model	Stand		Objective		Plate		Illuminator		
6A	Pole		Rotated 1x- 2x 1x-3x 2x-4x		White/Black		No		
6A-W	Arm Curved				White/Black		No		
6C	Pole				Frosted Glass		Incident & Transmitted		
6C-W	Arm Curved		Frosted Glass		Incident & Transmitted				
Technical Data									
Objective	Eyepiece								Working Distance (mm)
	5X		10X		15X		20X		
	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	
1X	5X	20.0	10X	20.0	15X	15.0	20X	10.0	76
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5	
1X	5X	20.0	10X	20.0	15X	15.0	20X	10.0	76
3X	15X	7.0	30X	6.7	45X	5.0	60X	4.3	
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5	76
4X	20X	5.0	40X	5.0	60X	3.7	80X	3.2	



A22.1209-6A

A22.1209-6AW



A22.1209-6C

A22.1209-6CW





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China