



Jewelry Cordless LED Stereo Optical Microscope 10X / 15X A23.0903-BL3

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0903-BL3
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 6.7x - 45x
- Applications: Jewelry, Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes, stereo dissecting microscope**



Product Description

Jewelry Cordless Stereo Optical Microscope 6.7x - 45x A23.0903-BL3

Description:

- 0.67~4.5x Zoom Lens, with WF10x/22mm Eyepiece,
- Total Optical Magnification 6.7~45x
- Big Fan Pole Stand, Up & Bottom LED Light
- Fully Coated Optical System, Sharp & Clear Image

Specifications:

A23.0903-BL3 Zoom Stereo Microscope, 0.67x-4.5x Specification		
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	
Zoom Range	0.67~4.5x	
Zoom Ratio	1:6.5	
Eyepiece	WF10x/22mm, Both with Eye-cup	
Magnification	6.7~45x, Up to 2~225x with Auxiliary Lens	
Working Distance	100mm, With Standard 1x Lens. Up to 26~287mm With Auxiliary Lens	
Stand	Big Fan Pole Stand BL3, With Focusing Knob, Focusing Range 50mm	
Illumination	Up & Bottom LED Light	
Power	100~240V Wide Range Voltage	
A23.0903-BL3 Zoom Stereo Microscope, 0.67x-4.5x Optional Accessories		
Eyepiece	WF15x/15mm	A51.0923-15
	WF20x/10mm	A51.0923-20
	WF25x/9mm	A51.0923-25
Auxiliary Lens	0.3x, W.D. 287mm	A52.0928-03
	0.4x, W.D. 217mm	A52.0928-04
	0.5x, W.D. 177mm	A52.0928-05
	0.75x, W.D. 117mm	A52.0928-75
	1.5x, W.D. 47mm	A52.0928-15
	2.0x, W.D. 26mm	A52.0928-20



A23.0903-BL3 Binocular Stereo Microscope



A23.0903-BL3 Binocular Stereo Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China