



WF10x - 20mm Clinic Stereo Optical Microscope Digital Stereo Microscopes A23.0905-BL8

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0905-BL8
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 8-35x
- Applications: Industry, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**



Product Description

Description:

A23.0905-BL8 Zoom Stereo Microscope, 0.8x-3.5x, Designed for common industry inspection and biological observation, greenough optical, magnification 8~35x. Zoom stereo len's magnification 0.8~3.5x, wide field, high eyepoint WF10x/20mm eyepiece.

--0.8~3.5x Zoom Lens, with WF10x/20mm Eyepiece,

--Total Optical Magnification 8~35x

--Pole Stand, Up & Bottom Light Source

Clinic Stereo Optical Microscope Digital Stereo Microscopes A23.0905-BL8

Feature :

A23.0905-BL8 Zoom Stereo Microscope, 0.8x-3.5x, Designed for common industry inspection and biological observation, greenough optical, magnification 8~35x. Zoom stereo len's magnification 0.8~3.5x, wide field, high eyepoint WF10x/20mm eyepiece.

--0.8~3.5x Zoom Lens, with WF10x/20mm Eyepiece,

--Total Optical Magnification 8~35x

--Pole Stand, Up & Bottom Light Source

Specifications:

A23.0905-BL8 Zoom Stereo Microscope, 0.8x-3.5x Specification		
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	
Zoom Range	0.8~3.5x	
Zoom Ratio	1:4.4	
Eyepiece	WF10x/20mm	
Magnification	8~35x	
Working Distance	97mm, With Standard 1x Lens	
Stand	Pole Stand BL8, With Focusing Knob, Stage Height 275*180*47mm	
Illumination	Up & Bottom Light Source	
Power	100~240V Wide Range Voltage	
A23.0905-BL8 Zoom Stereo Microscope, 0.8x-3.5x Optional Accessories		
Eyepiece	WF15x/15mm	A51.0923-15
	WF20x/10mm	A51.0923-20



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China