



Trinocular / Binocular Head Stereo Optical Microscope Rohs A23.0906-B4

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0906-B4
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7-45x
- Applications: Industry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**



for more products please visit us on cnoec.com

Product Description

Trinocular / Binocular Stereo Optical Microscope Rohs A23.0906-B4

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Pole Stand, No Light Source

Specifications:

A23.0906-B4 Zoom Stereo Microscope, 0.7x-4.5x Specification		A23.0906-B4	A23.0906-B4T
Head	Binocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm	o	
	Trinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm		o
Zoom Range	0.7~4.5x		
Zoom Ratio	1:6.5		
Eyepiece	WF10x/20mm, Both with Eye-cup		
Magnification	7~45x, Up to 3.5~180x with Auxiliary Lens		
Working Distance	100mm, With Standard 1x Lens		
Field of View	1~95.2mm		
Stand	Pole Stand B4, With Focusing Knob		
Illumination	No Light Source		
A23.0906-B4 Zoom Stereo Microscope, 0.7x-4.5x Optional Accessories			
Eyepiece	WF15x/15mm	A51.0922-15	
	WF20x/10mm	A51.0922-20	
Auxiliary Lens	0.5x	A52.0927-05	
	0.75x	A52.0927-75	
	1.5x	A52.0927-15	
	2.0x	A52.0927-20	
LCD Monitor	Add A59.1130 LCD Camera, it will became A36.0901-B4 LCD Stereo Microscope		



A23.0906-B4T
Trinocular Stereo Microscope



A23.0906-B4 Binocular Stereo Microscope



A23.0906-B4T
Trinocular Stereo Microscope



A23.0906-B4 Binocular Stereo Microscope

www.cnoec.com

www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

