



8X - 56X Medical Stereo Optical Microscope , CE Rohs A23.0907-B4

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0907-B4
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Magnification: 8X - 56X
- Applications: Industry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

8X - 56X Medical Stereo Optical Microscope , CE Rohs A23.0907-B4

Description:

A23.0907-B4 Zoom Stereo Microscope, 0.8x-5.6x, Designed for common industry inspection and biological observation, greenough optical, magnification 8~56x. Zoom stereo len's magnification 0.8~5.6x, Plan 10x/22mm eyepiece.

--Parallel Optical System, With Two Independent Optial Path

--0.8~5.6x Zoom Lens, with PL10x/22mm Eyepiece,

--Total Optical Magnification 8~56x

--Pole Stand, No Light Source

--Fixed Position and Fixed Times Function


--With Simple Operation, Excellent Optical Performance

Specifications:

A23.0907-B4 Zoom Stereo Microscope, 0.8x-5.6xSpecification	
Optical System	Parallel Optical System, With Two Independent Optial Path
Head	Brinocular, 45° Inclined, 360° Rotatable, +/-5° Diopter Adjustable, Interpupillary Adjustment 54~75mm
Zoom Range	0.8~5.6x
Zoom Ratio	1:7
Eyepiece	PL10x/22mm
Magnification	8~56x
Working Distance	90mm, With Standard 1x Lens
Stand	Pole Stand B4, With Focusing Knob
Illumination	No Light Source




A23.0907-B4 Binocular Stereo Microscope



A23.0907-B4 Binocular Stereo Microscope System Size



A23.0907-B4 Binocular Stereo Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China