



## Inspection Gem Stereo Optical Microscope 10x - 40x A23.1201-EW

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1201-EW
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 10x - 40x
- Applications: Gem, Jewelry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,  
stereo dissecting microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Inspection Gem Stereo Optical Microscope 10x - 40x A23.1201-EW

#### Description:

- Sharp stereo erect image over a wide field view.
- 45 degree inclined and 360° rotational binocular head with locked-eyepiece. The left ocular-tube with diopter adjustment  $\pm 5$ dp, The interpupillary distance is between 55-75mm
- Zoom objective 1x-4x
- Electric apparatus has CE approval
- Power supply 110V/60Hz or 220V/50Hz
- Each set packed in styrofoam and in an inner box. Two sets in one carton
- A wide range of optional accessories.

#### Specifications:

A23.1201-EW Specification										
Model	Head	Eyepiece	Objective	Stand	Working Stage					
A23.1201-EW	Binocular	SWF10X	Zoom1x-4x	Arm curved stand	Frosted Glass					
A23.1201-EW-3	Trinocular	SWF10X	Zoom1x-4x	Arm curved stand	Frosted Glass					
Technical Data										
Eyepiece	Standard Objective 1-4X		Auxiliary Objective							
			0.5X		0.75X		1.5X		2.0X	
	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance
5X	20-5	91	40-11	170	25-6.5	100	13-3.5	42	10-2.5	28
10X	23-5.5	91	52-12	170	30-7	100	15-4.8	42	11-5.3	28
15X	15.5-4	91	36-8.5	170	21-5	100	10.5-2.5	42	8-2	28
20X	10.5-3	91	25-5.8	170	14-3.5	100	7-1.8	42	5.5-1.5	28
Optional Accessories										Item No.
Eyepiece	WF5X /22									A51.1201-05
	WF10X /23									A51.1201-10
	WF15X/13									A51.1201-15
	WF20X/10									A51.1201-20
Auxiliary Objective	0.5X									A52.1282-05
	0.75X									A52.1282-75
	1.5X									A52.1282-15
	2X									A52.1282-20
Figure Eyepiece	VD-TV									A51.1207-TV
	VD-PC									A51.1207-PC
Stage	Darkfield stage, $\Phi 95$ mm									A54.1202-1
Plate	White/black Working Plate $\Phi 95$ mm									A54.1203-2
Illumination	more choice pls click here:light source									
Packing	Wooden cabinet with lock and key									A50.1202
	Aluminium case									A50.1203



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China