



## Small Handheld Monocular Inverted Optical Microscope A14.1020

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.1020
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Inverted Optical Microscope
- Theory: Inverted Microscope
- Highlight: **inverted biological microscope,**  
**inverted light microscope**



#### More Images




## Product Description

### Description:

- Small hand-held microscope, simple, good performance, high-definition.
- Inverted-design is available to observe liquid samples.
- Compact design, owned power, easy to carry out to take wildlife observation and inspect protozoa in factory effluent and microbiology.
- Small Size for easy portable: 14\*4\*18.5CM

### Specifications:

A14.1020 Inverted Biological Microscope		
Eyepiece	10x/ Φ18mm	
Objective	Achromatic 4x/0.08, W.D.16mm	
	Achromatic 10x/0.25, W.D.3.6mm	
	Achromatic 40x/0.6, W.D.0.6mm	
Nosepiece	Triple	
Stage	Plain Stage with Slide Clips	
Range of Vertical Stroke	Up 1.6mm, down 1.0mm, total 2.6mm	
Illumination	Brightfield, LED Illumination, AA Battery * 2	
A14.1020 Inverted Biological Microscope, Optional Accessories		
Objective	Achromatic 20x/0.42, W.D.2.85mm	A52.1002-20
Phase Contrast	Phase Contrast 10x	A5C.1022-10
	Phase Contrast 20x	A5C.1022-20
	Phase Contrast 40x	A5C.1022-40
	Phase Contrast, LED Illumination, AA Battery *2	A5C.1022-LED
Box	Portable Carring Box	A50.1006





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China