



1000x White Inverted Metallurgical Optical Microscope Research A13.2601

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2601
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research, Metallurgical
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 100x - 1000x
- Theory: Metallurgical Optical Microscope
- Highlight: **reflected light microscope,**
bright field microscope



Product Description

1000x White Inverted Metallurgical Optical Microscope Research A13.2601

Description:

A13.2601 Metallurgical Microscope, new design with high quality optical system, supplying clear image, is good to be used in laboratory research work or school education field.

- Plan Metallurgical Objective
- LWD Plan Metallurgical Objective Optional
- Polarizing Set of Polarizer & Analyzer
- Mechanical, size 180x165mm, moving range 50x40mm, 3 Stage Plate Available

Specifications:

Specification	A13.2601 Metallurgical Microscope		
Head	Trinocular, inclined 45°, interpupillary distance 50-75mm, light distribution 20:80		
Eyepiece	WF10X/18mm		
	WF12.5X/14mm		
	WF10x/20mm with micrometer 0.1mm/Div, 1pc		
Objective	Plan Metallurgical Objective		
	10x/0.25, W.D.=8.8mm		
	20x/0.35(S), W.D.=1.59mm		
	40x/0.65(S), W.D.=0.73mm		
	100x/1.25(O,S), W.D.=0.69mm		
Nosepiece	Quadruple nosepiece		
Focusing	Coaxial coarse and fine focusing knob		
Stage	Mechanical, size 180x165mm, moving range 50x40mm		
	Stage plate 1 Φ10		
	Stage plate 2 Φ20		
Polarizing	Polarizer Slide		
	Analyser Slide		
Filter	Blue, Green, Grey, White		
Light source	6V 30W halogen lamp with brightness control, spare lamp 2 pcs, spare fuse 2 pcs		
	Input Voltage 100-240V		
Packing	1pc/carton, size 59x59x33cm, G/N weight 15/11 Kgs		
Optional Accessories			Item No.
Eyepiece	Huygens 5x/20mm		A51.2601-05
	WF10x/20mm		A51.2602-1020
	WF16x/13mm		A51.2602-1613A
	WF10x/18mm, with micrometer 0.1m/Div		A51.2605-10
Objective	Plan Metallurgical Objective	2.5x/0.07	A5M.2630-2.5
		4x/0.10, W.D=17.5mm	A5M.2630-4
	LWD Plan Metallurgical Objective	PLL40x/0.60, W.D.=3.64mm	A5M.2631-40P
		PLL50x/0.70, W.D.=3.01mm	A5M.2631-50P
		PLL80x/0.80, W.D.=1.08mm	A5M.2631-80P
Stage Plate	Stage plate 3 Φ40		A54.2603-40
Adapter	CCD Adapter 0.6x		A55.2601-M06
	CCD Adapter 1.0x		A55.2601-M10
Digital Camera Adapter	For OLYMPUS digital camera		A55.2620
	For CANON digital camera		A55.2621
Digital SLR Camera Adapter	For NIKON SLR digital camera		A55.2622-N
	For SONY SLR digital camera		A55.2622-S
	For Canon SLR digital camera		A55.2622-C



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China