



## 100X - 1000X Research Inverted Metallurgical Microscopes A13.2602

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2602
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Research, Metallurgical
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular
- Illumination: Halogen Lamp, LED
- Magnification: 100X - 1000X
- Theory: Metallurgical Optical Microscope
- Highlight: **reflected light microscope, halogen lamp microscope**



## Product Description

### 100X - 1000X Research Inverted Metallurgical Microscopes A13.2602

#### Description:

**A13.2602 Metallurgical Microscope**, new design with high quality optical system, supplying clear image, is good to be used in laboratory research work or school education field.

- Compact Metallurgical Microscope at competitive cost
- Plan Metallurgical Objective
- LWD Plan Metallurgical Objective Optional
- Mechanical Stage Size: 180x180mm, Travel Range: 50\*40mm, 3 Stage Plate Available

#### Specifications:

Specification		A13.2602 Metallurgical Microscope		A	B
Head	Monocular head Inclined 45°			●	
	Binocular, inclined 45°, interpupillary distance 50-75mm				●
Eyepiece	WF10X/18mm			●	●
Objective	Plan Metallurgical Objective				
	10x/0.25, W.D.=8.8mm			●	●
	40x/0.65(S), W.D.=0.73mm				
	100x/1.25(O,S), W.D.=0.69mm				
Stage	Mechanical Stage Size: 180x180mm, Travel Range: 50*40mm			●	●
	Stage plate 1 Φ10			●	●
	Stage plate 2 Φ20			●	●
Focusing	Coaxial coarse and fine focusing knob			●	●
Light Source	3W LED Lamp, Input voltage:100V 240V			●	●
Filter	Blue, Green, Grey			●	●
Packing	1pc/carton, size: 44x27x52cm, G/N weight: 9/7Kgs			●	●
Optional Accessories			Item No.		
Eyepiece	Huygens 5x/20mm			A51.2601-05	
	WF10x/20mm			A51.2602-1020	
	WF12.5x/14mm			A51.2602-125	
	WF16x/13mm			A51.2602-1613A	
	WF10x/18mm, with micrometer 0.1m/Div			A51.2605-10	
Objective	Plan Metallurgical Objective	2.5x/0.07		A5M.2630-2.5	
		4x/0.10,W.D=17.5mm		A5M.2630-4	
		20x/0.35,W.D=1.59mm		A5M.2630-20	
	LWD Plan Metallurgical Objective	PLL40x/0.60, W.D.=3.64mm		A5M.2631-40P	
		PLL50x/0.70, W.D.=3.01mm		A5M.2631-50P	
		PLL80x/0.80, W.D.=1.08mm		A5M.2631-80P	
Illumination	6V/20W Halogen Lamp			A56.2613-6V20W	
Stage Plate	Stage plate 3 Φ40			A54.2603-40	



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China