



LWD Infinity Plan Metallurgical Optical Microscope Trinocular A13.2603

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2603
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research, Metallurgical
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Magnification: 50x - 500x
- Theory: Metallurgical Optical Microscope
- Highlight: **reflected light microscope,
bright field microscope**



Product Description

LWD Infinity Plan Metallurgical Optical Microscope Trinocular A13.2603

Description:

A13.2603 Metallurgical Microscope

- Infinity Plan Optical System
- Infinity Plan & LWD Infinity Plan Metallurgical Objective (no cover glass)
- Reflecting light, 6V 20W & 6V30W halogen lamp,brightness control
- Five-hole filter unit (blue, green, amber, grey, empty)

Specifications:

Specification	A13.2603 Metallurgical Microscope		A	B	C	D
Head	Binocular, inclined 30°, interpupillary distance 50-75mm		●		●	
	Trinocular, inclined 30°, interpupillary distance 50-75mm, light split 80:20			●		●
Eyepiece	WF10X/18mm, Reticule Eyepiece WF 10×-18 mm (0.1mm)		●	●	●	●
Objective	Infinity Plan Metallurgical Objective (no cover glass)	5x/0.11,W.D.=12.10	●	●	●	●
		10x/0.25,W.D.=4.75	●	●	●	●
	LWD Infinity Plan Metallurgical Objective (no cover glass)	20x/0.4,W.D.=8.35			●	●
		40x/0.65,W.D.=3.9	●	●		
		50x/0.7,W.D.=1.95			●	●
Stage	Mechanical Stage Size: 200×140mm, Travel Range: 76×52mm		●	●	●	●
Focusing	Coaxial coarse and fine focusing knob, coarse focusing range 25mm, fine focusing scale 2um, with tension adjustment		●	●	●	●
Nosepiece	Quadruple nosepiece		●	●	●	●
Light Source	Reflecting light, 6V 20W halogen lamp,brightness control		●	●		
	Reflecting light, 6V 30W halogen lamp,brightness control				●	●
Filter	Five-hole filter unit(blue, green, amber, grey, empty)		●	●	●	●
Packing	1pc/carton, size 56x31x41cm, G/N weight 9/7Kgs		●	●	●	●
Optional Accessories			Item No.			
Eyepiece	WF 5x/20mm		A51.2602-05			
	WF12.5x/14mm		A51.2602-125			
	WF16x/13mm		A51.2602-1613			
Objective	Infinity Plan Metallurgical Objective (no cover glass)	20x/0.4,W.D.=2.14	A5M.2632-20			
		40x/0.65,W.D.=3.9	A5M.2632-40			
	LWD Infinity Plan Metallurgical Objective (no cover glass)		80x/0.8,W.D.=0.85 A5M.2633-80			
Polarizing	Analyzer/ Polarizer		A5P.2603			
C-Mount	1x		A55.2601-MT10			
	0.47x, Focus Adjustable		A5P.2601-MT47			
	0.5x, Focus Adjustable		A5P.2603-MT50			



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China