



Trinocular Compensation Free Fluorescence Microscope Learning A16.0902

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.0902
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification


- Application: Learning|Laboratory
- Certification: CE|Rohs
- Drawtube: Binocular|trinocular
- Illumination: Halogen Lamp
- Magnification: 100X - 1000X
- Theory: Fluorescence Microscope
- Highlight: **epi-fluorescence microscope,**
inverted fluorescence microscope







Product Description

Specifications:

| A16.0902 Specification | | A | B |
|------------------------------------|---|-----------------|---|
| Head | Compensation free binocular head, 30 degree inclined (infinity) | o | - |
| | Compensation free trinocular head, 30 degree inclined (infinity) | - | o |
| Eyepiece | PL10X/22mm high eye-point plan eyepieceeach with eye-cap | | |
| Infinity Plan Achromatic Objective | Infinity plan achromatic objective 10X | | |
| | Infinity plan achromatic objective 210X | | |
| | Infinity plan achromatic objective 40X | | |
| | Infinity plan achromatic objective 100X | | |
| Nosepiece | Reversed quintuple nosepiece | | |
| Focus | Low adjustment coaxial focus system with limit stopper, moving range: 25mm, precision: 0.002mm | | |
| Stage | 175*145mm mechanical stage,adjusting for right hand | | |
| Condenser | N.A.1.2/0.22 swing achromatic condenser | | |
| Illuminator | 90V~240V/ 6V 30W halogen bulb,brightness adjustable | | |
| Fluorescent Kits | Power box: 100V-240VAC with Indicator Display | | |
| | 100WHBO Ultra Hi-voltage Spherical Mercury Lamp OSRAM | | |
| | B(Blue)excitation filter system | | |
| | G(Green)excitation filter system | | |
| Optional Accessories | | | |
| Eyepiece | High eye-point PL10x/22mm with micrometer, diopter adjustable | A51.0905-1022R | |
| | High eye-point PL15x/16mm | A51.0903-1516 | |
| Objectvie | Infinity plan objective 4x | A5F.0931-4 | |
| Stage | 175*145mm mechanical stage,adjusting for left hand | A54.0941-10 | |
| Color Filter | Yellow (Dia:45mm) | A56.0935-XY-FCY | |
| | Green(Dia:45mm) | A56.0935-XY-FCG | |
| | Blue(Dia:45mm) | A56.0935-XY-FCB | |
| | Neutral color filter (Dia:45mm) | A56.0935-XY-FCF | |
| CCD adapter | 1x C-Mount,Focus Adjustable | A55.0926-10 | |
| | 0.35x, For 1/3 C-Mount,Focus Adjustable | A55.0926-35 | |
| | 0.50x, For 1/2 C-Mount,Focus Adjustable | A55.0926-50 | |
| | 0.67x, For 2/3 C-Mount,Focus Adjustable | A55.0926-67 | |







Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China