



## Petrology / Geology Polarizing Light Microscope Binocular A15.1014

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1014
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Research, Laboratory, Metal, Silicon Wafers
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 40X - 600X
- Theory: Polarizing Light Microscope
- Highlight: **student compound microscope, compound light microscope**



#### More Images



## Product Description

### Specifications:

A15.1014 Polarizing Microscope Specification		M	B	T
Viewing head	Monocular Head, inclined at 45°, 360° rotatable	●	-	-
	Compensation Free binocular head at 30°, 360° rotatable		●	-
	Compensation Free Trinocular head at 30°, 360° rotatable	-		●
Eyepiece	Eyepiece WF10x /18 with scale of crosshair	●	●	●
	WF10X/18	-	●	●
Objective	Strain free Achromatic 4X, 10X, 40X, 60X	●	●	●
Nosepiece	Backward quadruple nosepiece	●	●	●
Analyzer	Rotatable analyzer with gradation 0°-90°	●	●	●
Bertrand Lens	Sliding in/ out of Optical path	●	●	●
Optical Compensation	$\lambda$ Slip (first class red),	●	●	●
	1/4 $\lambda$ Slip quartz wedge	●	●	●
Polarizer	Sliding in/out of Optical Path , Located on the Top of Collector	●	●	●
Condenser	N.A.1.25 Abbe Condenser with Iris Diaphragm & Filter	●	●	●
Stage	Diameter $\Phi$ 145mm, 360° rotatable and graduated in 1° increments, minimum resolution 6' when using vernier scale, center adjustable	●	●	●
Focusing	Coaxial Coarse and Fine Focusing System, Range 28mm, Division 0.002mm	●	●	●
Illumination	Halogen Lamp 220V or 110V/6V20W, Brightness Adjustable	●	●	●



[www.cnoec.com](http://www.cnoec.com)





A15.1014-B



A15.1014-M

A15.1014-T



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China