



40x - 1000x Polarizing Light Microscope Halogen Lamp A15.1018

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1018
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research, Laboratory, Metal
- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Theory: Polarizing Light Microscope
- Highlight: **student compound microscope ,
compound light microscope**



Product Description

Specifications:

A15.1018 Polarizing Microscope Specification		
Optical System	Infinite optical system	
Head	Seidentopf trinocular head inclined at 30° , Interpupillary 48-75mm	
Eyepiece	Plan eyepiece EW10x/22	
	Plan eyepiece with scale of crosshair, Plan eyepiece with cross, Plan eyepiece with gridding	
Objective	Strainfree achromatic objective 4x/0.10/∞/- WD 17.3mm	
	Strainfree achromatic objective 10x/0.25/∞/- WD 10.0mm	
	Strainfree achromatic objective 40x/0.65/∞/- WD 0.54mm	
	Strainfree achromatic objective 100x/1.25/∞/- WD 0.13mm	
Nosepiece	Backward quintuple nosepiece, center adjustable	
Analyzer	360° dial rotation Minimum scale reading: 0.1°(Vernier scale)	
Conoscopic Observation	Switch between orthoscopic and conoscopic observation, Bertrand lens position adjustable	
Optical Compensator	--λSlip(first class red)	
	--1/4λSlip	
	--Quartz wedge	
Polarizing	Dia Φ170mm, center adjustment, 360° scale, minimum division 1° ,	
Rotatable Stage	minimum reading 6' by means of vernier scale, 45° click stop knob	
Condenser	Strain-free swing condenser NA.0.9/0.25	
Focusing	Coaxial coarse and fine adjustment , range 32mm , fine division 0.001 mm, coarse stroke 37.7mm per rotation, fine stroke 0.1mm per rotation	
Polarizing	360° rotatable 0° position adjustable	
Filter	Φ45mm blue filter	
Illumination	Halogen Lamp 24V/100W, brightness adjustable	
Optional Accessories		Item No.
Objective	Strain-free Plan Achromatic objective 20x/0.40/∞/- WD 5.10mm	A5P.1030-20
Attachable Stage	Moving range 30mmx40mm	A54.1002-T
Filter	Φ45mm ND6 Neutral density filter	A5F.1011-6
	Φ45mm ND25 Neutral density filter	A5F.1011-25
Attachment	Video attachment	A55.1002



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

