



Lab Dual Head Multi Viewing Microscope 40x - 1000x A17.1013-A

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1013-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Lab|educational|school
- Certification: CE|Rohs
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Drawtube: Binocular
- Highlight: **dual head microscope, Multi View Microscope**



Product Description

Description:

Multi Viewing Microscope, through a multiple optical path, allows 2~6 people observation at the same time. With the features of high quality optical system, effective Illumination, LED pointer, the main operator can easily show observation and process to colleagues or students. Coherence, it is widely used in Clinic, Scientific Research and Teaching Demonstration

Specifications:

A17.1013-A Multi-Viewing Microscope, 2 Position	
Head	2 pcs Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm
Eyepiece	WF10X, 4 pcs
Objective	Achromatic 4x, 10x, 40x(s), 100x(s, oil)
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double Layer Mechanical Stage 140mm x 140mm, Travel Range 75x50mm
Condenser	Abbe N.A.=1.25 with Iris Diaphragm & Filter
Focusing	Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.2mm Per Rotation, Moving Range 28mm
Light Source	Base with Built-in Halogen Lamp 6V/20W , with Continue Intensity Control, 110V/220V

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China