



## 1000x Multi Viewing Microscope Penta Head Microscopes A17.1013-C

### Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1013-C
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Educational, School
- Certification: CE, Rohs
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Drawtube: Binocular
- Highlight: dual head microscope, Multi View Microscope



#### More Images



## Product Description

### Specifications:

|                                                 |                                                                                                                                                   |
|-------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| A17.1013-C Multi-Viewing Microscope, 5 Position |                                                                                                                                                   |
| Head                                            | 5 pcs Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm                                                                        |
| Eyepiece                                        | WF10X, 10 pcs                                                                                                                                     |
| Objective                                       | Achromatic 4x, 10x, 40x(s), 100x(s, oil)                                                                                                          |
| Nosepiece                                       | Quadruple Revolving Nosepiece                                                                                                                     |
| Stage                                           | Double Layer Mechanical Stage 140mm x 140mm, Travel Range 75x50mm                                                                                 |
| Condenser                                       | Abbe N.A.=1.25 with Iris Diaphragm & Filter                                                                                                       |
| Focusing                                        | Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.2mm Per Rotation, Moving Range 28mm |
| Light Source                                    | Base with Built-in Halogen Lamp 24V/100W , with Continue Intensity Control, 110V/220V                                                             |
| Indicator                                       | Demo. Indicator with Illuminator of 6V/3W                                                                                                         |





A17.1013 Multi-Viewing Microscope



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China