



2 Position , Infinity Plan Multi Viewing Microscope With Kohler Illumination A17.1025

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1025
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Lab, Educational
- Certification: CE, Rohs
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Highlight: dual head microscope, Multi View Microscope



Product Description

2 Position , Infinity Plan Multi Viewing Microscope With Kohler Illumination A17.1025

Specifications:

| | | |
|---|---|-------------|
| A17.1025 Multi-Viewing Microscope, 2 Position | | |
| Optical System | Infinite Optical System | |
| Head | 1 pcs Trinocular Head, 30°Inclined, 360° Rotatable, Interpupillary Distance 48-75mm | |
| | 1 pcs Brinocular Head, 30°Inclined, 360° Rotatable, Interpupillary Distance 48-75mm | |
| Eyepiece | WF10X, 4 pcs | |
| Objective | Infinity Plan Achromatic 4x, 10x, 40x(s), 100x(s, oil) | |
| Nosepiece | Reversed Quadruple Revolving Nosepiece | |
| Stage | Double Layers Mechanical Stage 185x142mm, Moving Range 75x55mm | |
| Focusing | Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm | |
| Condenser | Swing Abbe Condenser N.A.0.9/0.25 | |
| Light Source | External Kohler Illumination, Halogen Lamp 12V50W | |
| Pointer | Green LED Pointer, Brightness Adjustable | |
| A17.1026-A Multi-Viewing Microscope, Optional Accessories | | |
| Adapter | 1.0x C-Mount | A55.1002-10 |
| | 0.5x C-Mount | A55.1002-05 |



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China