



Educational Multi View Microscope Achromatic Microscopes CE A17.1026-C

Our Product Introduction

for more products please visit us on cnoec.com


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1026-C
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research|Lab|educational
- Certification: CE
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Highlight: Multi View Microscope, Dual View Microscope



More Images




Product Description

Our Product Introduction

for more products please visit us on cnoec.com

Specifications:

A17.1026-C Multi-Viewing Microscope, 5 Position		
Optical System	Infinite Optical System	
Head	1 pcs Trinocular Head, 30°Inclined, 360° Rotatable, Interpupillary Distance 48-75mm 4 pcs Brinocular Head, 30°Inclined, 360° Rotatable, Interpupillary Distance 48-75mm	
Eyepiece	WF10X, 10 pcs	
Objective	Infinity Plan Achromatic 4x, 10x, 40x(s), 100x(s, oil)	
Nosepiece	Reversed Quadruple Revolving Nosepiece	
Stage	Double Layers Mechanical Stage 185x142mm, Moving Range 75x55mm	
Focusing	Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm	
Condenser	Swing Abbe Condenser N.A.0.9/0.25	
Light Source	External Kohler Illumination, Halogen Lamp 12V50W	
Pointer	Green LED Pointer, Brightness Adjustable	
A17.1026-C Multi-Viewing Microscope, Optional Accessories		
Adapter	1.0x C-Mount	A55.1002-10
	0.5x C-Mount	A55.1002-05





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China