



## Lab Research Dual Mult Viewing Microscope Halogen Lamp 3 Position A17.1102-A

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1102-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Application: Research|Lab|educational
- Certification: CE
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Highlight: dual head microscope, Dual View Microscope



## Product Description

### Specifications:

A17.1102-A Multi-Viewing Microscope, 3 Position	
Head	3 pcs Compensation Free Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm
Eyepiece	WF10X/22mm, 6 pcs
Objective	Infinity Achromatic 4x, 10x, 40x(s), 100x(s, oil)
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double Layer Mechanical Stage 180mm x 150mm, Travel Range 75x50mm
Condenser	Abbe N.A.=1.25 with Iris Diaphragm & Filter
Focusing	Coaxial Coarse & Fine Focusing Ddjustmen Fine Fousing Scale Value 0.002mm
Light Source	Halogen Bulb 12V30W. AC85V-230V Adjustable Brightness, High Light Kohler Illumination

### Description:

Multi Viewing Microscope, through a multiple optical path, allows 2~6 people observation at the same time. With the features of high quality optical system, effective Illumination, LED pointer, the main operator can easily show observation and process to colleagues or students. Coherence, it is widely used in Clinic, Scientific Research and Teaching Demonstration

### About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China