



## Black Inclined Binocular Jewelry Microscope 20X - 40X A24.1201

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 20X - 40X
- Theory: Jewelry Microscope
- Highlight: **laboratory binocular microscope, gemological microscope**



## Product Description

### Description:

Jewelry Microscope, is special designed microscope on the base of stereo microscope, with jewelry clamp, dark field light source, and more illumination ways, for jewelry and gem inspection and observation.

### Specifications:

A24.1201 Jewery Microscope		
Head	Inclined Binocular Head,360° Rotatable	
Eyepiece	WF10x	
Objective	2x4x	
Incident Light	Halogen Lamp 12V/10W, Brightness Adjustable	
Transmitted Light	Halogen Lamp 12V/10W, Brightness Adjustable	
Working Stage	Jewel Tweezer	
	Dark Field Working Stage, With Diaphram	
A24.1201 Jewery Microscope , Optional Accessories		
Eyepiece	WF5x	E51.1221-05
	WF15x	E51.1221-15
	WF20x	E51.1221-20



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China