



Black Inclined Binocular Jewelry Microscope 20X - 40X A24.1201

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1201
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 20X - 40X
- Theory: Jewelry Microscope
- Highlight: **laboratory binocular microscope,
gemological microscope**



Product Description

Description:

Jewelry Microscope, is special designed microscope on the base of stereo microscope, with jewelry clamp, dark field light source, and more illumination ways, for jewelry and gem inspection and observation.

Specifications:

A24.1201 Jewery Microscope		
Head	Inclined Binocular Head,360° Rotatable	
Eyepiece	WF10x	
Objective	2x4x	
Incident Light	Halogen Lamp 12V/10W, Brightness Adjustable	
Transmitted Light	Halogen Lamp 12V/10W, Brightness Adjustable	
Working Stage	Jewel Tweezer Dark Field Working Stage, With Diaphram	
A24.1201 Jewery Microscope . Optional Accessories		
Eyepiece	WF5x	E51.1221-05
	WF15x	E51.1221-15
	WF20x	E51.1221-20



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China