



20W 50W LED Cold Light Source Microscope Accessories Rohs A56.0602

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.0602
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Voltage: AC 100V~240V 50Hz~60Hz
- Power: 20W,50W
- Color T: 6400K
- Highlight: compound microscope parts,
parts of the microscope



Product Description

Specifications:

	A56.0602 LED Cold Light Source Specification	A	B
Voltage	AC 100V~240V 50Hz~60Hz		
Power	20W	•	
	50W		•
Color T	6400K		
Luminous flux	>280 lm	•	
	>500 lm		•
Dimension	210×450×80(mm)		
N. W.	1.6 Kg		
Lamp Life	50000h		
Cooling	convection cooling		
Brightness Control	Continues		
Dual Pipe Light Guide	Fiber Optic Goose Neck Dual Pipe With Focusing Lens, Fiber Dia. 1/4" (6mm), 500mm In Length, 5/8" (16mm) Standard Adapter For Attaching The Cold Light Source		



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China