



## 2700K / 6000K microscope Fluorescent ring light Microscope Accessories A56.1203

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1203
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Voltage: 110V/220V
- Power: 10W
- Color Temperature: 2700K/6000K
- Highlight: **compound microscope parts,  
parts of the microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Specifications:

Model	A56.1203-1	A56.1203-2	A56.1203-2T
Type	No Light Adjustment		Light Adjustment
Bulb	A56.1203-L1 2 Pins, 10W	A56.1203-L2 4 Pins, 10W	
Voltage	110V/220V		
Power	10W		
Color Temperature	2700K/6000K		



A56.1203-1 Fluorescent Ring Light, 10W,  
Short Type with Aluminum Rotate Ring



A56.1203-L1 Fluorescent Lamp, 10W, 2 pin



A56.1203-2 Fluorescent Ring Light, 10W



A56.1202-L2 Fluorescent Lamp, 10W, 4 pin



A56.1203-2T Fluorescent Ring Light, 10W, Brightness Adjustable



A56.1202-L2 Fluorescent Lamp, 10W, 4 pin



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China