



24V 150W Optical Fiber Microscope Light Source A56.1223

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1223
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE
- Output: USB3.0
- Product Name: Microscope Accessories
- Power: 24V 150W
- Model Number: A56.1223
- Highlight: compound microscope parts,
parts of the microscope



More Images



Product Description

Description:

- Modular design for easy replacement of bulb, socket and computer board.
- Quick release of fiber optic cables for multiple light application requirements.
- Built-in aluminum alloy guard plate to prevent excessive heat from hurting the user.

Specifications:

A56.1223 Optical Fiber Cold Light Source	
A56.1223	Optical Fiber Cold Light Source Main Body 24V150W
A56.1223-1	Single Pipe Light Guide Length 550mm, Optical Fiber Dia. 5.6mm
A56.1223-2	Dual Pipe Light Guide Length 550mm, Optical Fiber Dia. 5.6mm
A56.1223-3	Ring Light Guide Ring Dia. 55.5mm, Length 800mm, Optical Fiber Dia. 8mm
A56.1223-3	Ring Light Guide Ring Dia. 55.5mm, Length 900mm, Optical Fiber Dia. 10mm



A56.1223 Fiber Optical Cold Light Source



A56.1223-1 Single Pipe Light Guide



A56.1223-2 Dual Pipe Light Guide



A56.1223-3, A56.1223-4 Ring Light Guide



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

