



Mono Zoom Video Microscope led Lamp Microscope Accessories A56.2102

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.2102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Inner Diameter: 1.6", 2.4"
- LED Power: 6W, 8W
- Number Of LEDs: 8
- Highlight: microscope accessory, parts of the microscope



More Images



for more products please visit us on cnoec.com

Product Description

Our Product Introduction

Description:

UV Free
 8-points high quality LED light
 1 watt per bulb
 Adjustable/dimming control
 ESD safe
 RoHS Compliant
 Build in optimum angle for even distribution of light
 Optional polarizer available to eliminate glare
 Long life LED bulbs

Specifications:

	A56.2102-1	A56.2102-2
Used For	Mono Zoom Video Microscope	Normal Stereo Microscope
Inner Diameter	1.6"	2.4"
Luminosity Brightness	1.9kfc	2.2kfc
LED Power	6W	8W
Color temperature	6000k~7000k (cool white light)	
Power input	96-265VAC	
Power output	30VDC	
Life	50,000h	
Number of LEDs	8	
Light Intensity Control	Yes	
Warranty	1 Year	
	Optional Accessory	
Polarizing Kits		Item No. A56.2102-PL



A56.2102-1 LED Lamp
for Mono Zoom Video Microscope

A56.2102-2 LED Lamp
for Normal Stereo Microscope



A56.2102-PL Polarizing Kit


for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China