



## Mono Zoom Video Microscope led Lamp Microscope Accessories A56.2102

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.2102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Inner Diameter: 1.6", 2.4"
- LED Power: 6W, 8W
- Number Of LEDs: 8
- Highlight: microscope accessory, parts of the microscope



#### More Images



## Product Description

### Description:

UV Free  
 8-points high quality LED light  
 1 watt per bulb  
 Adjustable/dimming control  
 ESD safe  
 RoHS Compliant  
 Build in optimum angle for even distribution of light  
 Optional polarizer available to eliminate glare  
 Long life LED bulbs

### Specifications:

	A56.2102-1	A56.2102-2
Used For	Mono Zoom Video Microscope	Normal Stereo Microscope
Inner Diameter	1.6"	2.4"
Luminosity	1.9kfc	2.2kfc
Brightness		
LED Power	6W	8W
Color temperature	6000k~7000k (cool white light)	
Power input	96-265VAC	
Power output	30VDC	
Life	50,000h	
Number of LEDs	8	
Light Intensity Control	Yes	
Warranty	1 Year	
Optional Accessory		Item No.
Polarizing Kits		A56.2102-PL



A56.2102-PL Polarizing Kit



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China