



OPTO-EDU A12.2601-A Laboratory Biological Microscope Binocular Achromatic

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.2601-A
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Eyepiece: WF10x/18mm
- Nosepiece: Quadruple, Backward
- Objective: Achromatic Objective
- Working Stage: Mechanical Stage Size: 145 Mm×140mm
- Focusing: Coaxial Coarse & Fine Focusing
- Condenser: Abbe Condenser N.A.1.25, Iris Diaphragm
- Illumination: 3W-LED Illumination Systems
- Highlight: **digital trinocular microscope, infinity plan microscope**



More Images



Product Description

Seidentopf Head

Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, High Eyepoint Eyepiece, Which Is Convenient For People With Different Sitting Height; The Integrated Design Of The Eyepiece Lens Barrel, Avoids The Loss And Contamination Of The Eyepiece And Enhances The Integrated Performance Of The Microscope



**Seidentopf Head Inclined 30°
Rotatable 360° , Interpupillary Distance 48-75mm**

Quadruple, Backward

4 Sets Objective Available: Achromatic, Plan, Infinity E-Plan, Infinity Plan.

Double Slide Holder Working Stage

Mechanical Stage Size: 145 mm×140mm, Travel: 76mmX52mm, Scale: 0.1mm



**Mechanical Stage Size 145×140mm
Travel Range 76×52mm**

Condenser & Light Source

Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder
3W-LED Illumination Systems



Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down



3W-LED Illumination Systems

Items	A12.2601 Biological Microscope Specification	-A	-AT	-B	-BT	-C	-CT	-D	-DT	Cata. No.
Head	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm	●	○	●	○	●	○	●	○	A53.2621-B
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20:E80	○	●	○	●	○	●	○	●	A53.2621-T

[illegible]



Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China