



## 200x Portable Usb Handheld Digital Microscope For Home , Rohs CE A34.4501

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Data Transfer: USB Connection Version 2.0
- Model Number: A34.4501
- Application: Research, Learning, Home
- Certification: CE, Rohs
- Magnification: 200x
- Theory: Handheld Digital Microscope
- Highlight: **Portable Digital Microscope,  
Hand held Digital Microscope**



#### More Images



## Product Description

### 20x Portable Usb Handheld Digital Microscope For Home , Rohs CE A34.4501

#### Description:

- Magnification 10x,60x,200x Changable!
- 1.3M Digital Microscope
- Up & Bottom Two Light Source Intergrated
- To Take Pictures and Record Videos
- With Processing Software Ulead Videostudio, Support 7 Langaues

#### Specifications:

|                                                                           |                                                           |
|---------------------------------------------------------------------------|-----------------------------------------------------------|
| A34.4501 USB Digital Microscope, 200x                                     |                                                           |
| Magnification                                                             | 10x, 60x, 200x Changable                                  |
| Data Transfer                                                             | USB Connection Version 2.0                                |
| Power Supply                                                              | USB Connection Version 2.0 (5v)                           |
| Image Resolution                                                          | RGB 1280x1024 pixels, 1.3M                                |
| Active Image Size                                                         | RGB 640x480 Pixels With Connection to USB2.0 Port         |
| Light Source                                                              | Up & Bottom Two LED Light Source Intergrated              |
| Software                                                                  | Ulead Videostudio 7SE, To Take Pictures and Record Videos |
| System Require                                                            | Windows XP, Windows 7 32 Bits                             |
| MOQ                                                                       | 10 Sets                                                   |
| A34.4501 USB Digital Microscope, Standard Outfit Box Packing List         |                                                           |
| USB Microscope                                                            |                                                           |
| USB Cable                                                                 |                                                           |
| Picture Editing Software Disc Ulead Videostudio 7SE(Not Shown in Picture) |                                                           |
| Box of 12 Sildes (5 Pre-made Dissect Slides, 7 Blank Slides)              |                                                           |
| 4 Sample Flasks (Adhesive, Eosin Dye, Sea Salt, Brine Shrimp Eggs)        |                                                           |
| Storage Container                                                         |                                                           |
| Culture Dish                                                              |                                                           |
| Plastic Tweezers                                                          |                                                           |
| Warranty Card                                                             |                                                           |



[www.cnoec.com](http://www.cnoec.com)

[www.cnoec.com](http://www.cnoec.com)



[www.cnoec.com](http://www.cnoec.com)





**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China