



Laboratory Binocular Compound Optical Microscope Phase Contrast A12.0205

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0205
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope, optical compound microscope**

Product Description

Laboratory Binocular Compound Optical Microscope Phase Contrast A12.0205


New Product in 2011!

A12.0205 biological microscopes is equipped with infinity plan achromatic objectives and wide field eyepieces. It have beautiful sculpt and provide clear image, wide field, convenience operation. It can be used for biology medicine agriculture and industry area widely. It is the ideal instrument in medical treatment, teaching demonstration and scientific research.

Specifications:

Specification	A12.0205	
Head	Trinocular, inclination of 30°	
Eyepiece	WF10X(Φ22mm)	
Infinity Plan Achromatic Objective	PL 4X/0.10 (Work distance):19.8 mm	
	PL 10X/0.25 (Work distance):5.0 mm	
	PL 40X/0.65 (Spring) (Work distance):0.66 mm	
	PL 100X/1.25 (Spring, oil) (Work distance): 0.36mm	
Focus System	Coaxial coarse/fine focus system, Minimum division of fine focusing: 2um With tensional adjustable and up stop,	
Nosepiece	Quadruple (backward ball bearing inner locating)	
Stage	Double layer mechanical, size 193x155mm, moving range 75x50mm	
Transmitted Illumination System	Abbe condenser NA.1.25 Rack & pinion adjustable	
	Blue filter and Ground glass	
	LED, High Brightness, Light Adjustable	
Optional Accessories		Item No.
Eyepiece	Dividing eyepiece(field number: μ22mm) , 0.1mm/Div	A51.0205-G10B
Infinity Plan Achromatic Objective	PL 20X/0.40 (Work distance:8.8 mm	A52.0206-20
	PL 60X/0.85 (Work distance):0.25mm	A52.0206-60
Filter	Green filter	A56.0209-02
	Yellow filter	A56.0209-03
Phase Contrast Unit	Centering telescope	A5C.0212-1
	Objective, Phase contrast plan 10x/0.25 PHP	A5C.0230-P10
	Objective, Phase contrast plan 20x/0.40 PHP(Spring)	A5C.0230-P20
	Objective, Phase contrast plan 40x/0.65 PHP(Spring)	A5C.0230-P40
	Objective, Phase contrast plan 100x/1.25 PHP (Spring,oil)	A5C.0230-P100
CCD adapter	0.4x	A55.0202-01
	0.5x	A55.0202-04
	1x	A55.0202-02
	0.5x with dividing 0.1mm/Div	A55.0202-03
Digital Camera Adapter	For Canon EF ditigal SLR camera adapter	A55.0204-02
Polarizing Attachment	Polarizer 360° rotatable	A5P.0210-1
	Analyzer	A5P.0211-1





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China