



40X - 1000X Ergonomic Binocular Head LED Light Source Microscope , A12.1029

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1029
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp, LED
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope ,
optical compound microscope**



Product Description

School Binocular LED Biological Microscope 40X - 1000X , A12.1029

Description:

Compound Optical Microscope, is a type of microscope which uses visible light and a system of lenses to magnify images of small samples. We have a wide range for usage in laboratory, university, institute, or schools, for used extensively in microelectronics, nanophysics, biotechnology, pharmaceuticals research, mineralogy and microbiology.

Specifications:

Specification A12.1029		A	B
Head	Compensation Free binocular head, 30° inclined, Interpupillary 48-75mm	●	
	Compensation Free trinocular head, 30° inclined, Interpupillary 48-75mm		●
Eyepiece	Wide Field Eyepiece WF 10×/18		
Nosepiece	Quadruple Nosepiece		
Objective	Infinitive Semi-plan Achromatic Objectives 4X		
	Infinitive Semi-plan Achromatic Objectives 10X		
	Infinitive Semi-plan Achromatic Objectives 40X		
	Infinitive Semi-plan Achromatic Objectives 100X		
Stage	Double Layers Mechanical Stage 140mm×140mm/ 75mm×50mm		
Condenser	Sliding-in centerable condenser N.A.1.25		
Focusing	Coaxial Coarse & Fine Adjustment, Moving Range 20mm		
Illumination	S-LED Lamp, Brightness Adjustable		
Optional Accessories		Item No.	
Eyepiece	Extra Wide Field Eyepiece EW10x/20 with Diopter Adjustment	A51.1005-1020	
Objective	Infinitive plan Achromatic Objectives 4X	A52.1006-4	
	Infinitive plan Achromatic Objectives 10X	A52.1006-10	
	Infinitive plan Achromatic Objectives 40X	A52.1006-40	
	Infinitive plan Achromatic Objectives 100X	A52.1006-100	
Illumination	6V/ 20W Halogen Lamp, Brightness Adjustable	A56.1003-1	
Phase Contrast Kit		A5C.1001	
Epi-fluorescent Attachment	XY-2 EPI-fluorescent Attachment	A5F.1001	
	FL-LED EPI-fluorescent Attachment		
Dark Field Attachment		A5D.1001	





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China