



Hobby Biological Compound Optical Microscope With Quadruple Nosepiece A12.1032

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1032
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory, Hobby, Home
- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **optical compound microscope, wide field microscope**



Product Description

Hobby Biological Compound Optical Microscope With Halogen Lamp A12.1032

Specifications:

Specification A12.1032		A	B
Head	Seidentopf Binocular Head, Inclined at 30°		
Eyeiece	Wide Field Eyepiece WF 10x/18		
Nosepiece	Quadruple Nosepiece		
Objective	Achromatic Objectives 4X		
	Achromatic Objectives 10X		
	Achromatic Objectives 40X		
	Achromatic Objectives 100X		
Stage	Double Layers Mechanical Stage 140mm×140mm/ 75mm×50mm		
Condenser	N.A.1.25 Abbe Condenser with Iris Diaphragm and Filter		
Filter	Blue Filter		
	Green Filter		
Illumination	6V/ 20W Halogen Lamp, Brightness Adjustable		
Other	Back Cover	-	●
Optional Accessories		Item No.	
Eyeiece	Wide Field Eyepiece WF16x/ 11	A51.1002-16	
Objective	Plan Achromatic Objectives 4X	A52.1004-4	
	Plan Achromatic Objectives 10X	A52.1004-10	
	Plan Achromatic Objectives 40X	A52.1004-40	
	Plan Achromatic Objectives 100X	A52.1004-100	
Polarization Attachment	Simple Polarization Set	A5P.1001	
Fluorescent Attachment	Epi-fluorescent Attachment	A5F.1001	
Photo Attachment		A55.1001	
Video Adapter	C Mount 1x	A55.1002-10	
	C Mount 0.5x	A55.1002-05	
Phase Contrast Attachment	Simple Phase Contrast Kit	A5C.1001	
Dark Field Attachment	Sliding Phase Contrast Kit	A5C.1002-C3	
	Turret Phase Contrast kit	A5C.1003-B3	
Dark Field Attachment	Dark Field Condenser (Dry) N.A. 0.9	A5D.1001-BD	
	Dark Field Condenser (Immersion) N.A. 0.9	A5D.1001-BO	





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China