



Monocular Head Biological Microscope 400X For Primary School Student A11.1111

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1111
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 400X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope, binocular biological microscope**



Product Description

Monocular Biological Microscope 400X For Primary School Student A11.1111

Features:

- This microscope is suitable for primary school student use
- Mechanical Length 160mm
- Full Coated Optical System.

Specification:

A11.1111 Biological Microscope		-A	-B
Head	Monocular Head,45° inclined	●	●
	Monocular Head,45° inclined,360° Rotatable	-	●
Eyepiece	WF10X	●	●
Objective	Achromatic 4X	●	●
	Achromatic 10X	●	●
	Achromatic 40X, Spring	●	●
Nosepiece	Triple Revolving Nosepiece	●	●
Stage	Plain Stage with Specimen Clips, Size 90mmx90mm	●	●
Condenser	Single Lens with Disc Diaphragm	●	●
Focusing	Coarse Focusing	●	
	Coaxial Coarse & Fine Focusing Adjustment Fine Focusing Scale Value 0.002mm	-	●
Light Source	Incandescent Lamp,220V/20W or 110V/20W	●	●
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
Objective	Achromatic 20X		
	Achromatic 60X(S)		
Light Source	LED		



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China