



Lab Achromatic Biological Microscope Monocular Phase Contrast Microscopes A11.1132

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1132
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 400X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,**
biological compound microscope

Product Description

Lab Achromatic Biological Microscope Monocular For Student A11.1132

Features:

Full metal body, high stability
Cordless, compact size, easy to carry
Top and bottom LED light source
Competitive low price

Specification:

A11.1132 Economical Cordless Student Microscope Specification	
Head	Monocular Head
Eyepiece	WF10X
Objective	Achromatic 4X
	Achromatic 10X
	Achromatic 40X, Spring
Nosepiece	Triple Revolving Nosepiece
Stage	Plain Stage with Clips, size 90*90mm
Condenser	6 Hole Disc Diaphragm
Focusing	Coarse Focusing Knobs
Light Source	Top and Bottom LED light, powered by AA 1.5V Battery 3 pcs. (No battery included when ship by air)
Body	Main Body Size 10*15*27cm

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



A11-1132



Triple Nosepiece, Upper LED Light



8 Hole Disc Diaphragm



Coarse Focusing Knobs



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China