



Binocular Stereo Optical Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1201-EC2
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 10X-40X
- Applications: Industry, Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



for more products please visit us on cnoec.com

Industry Binocular Stereo Optical Microscope A23.1201-EC2

- Sharp stereo erect image over a wide field view.
- 45 degree inclined and 360o rotational binocular head with locked-ineyepiece.The left ocular-tube with diopter adjustment $\pm 5dp$,The interpupillary distance is between 55-75mm
- Zoom objective 1x-4x
- Electric apparatus has CE approval
- Power supply 110V/60Hz or 220V/50Hz
- Each set packed in styrofoam and in an inner box.Two sets in one carton
- A wide range of optional accessories.

A23.1201-EC2 Specification											
Model	Head	Eyepiece		Objective		Stand		Working Stage			
A23.1201-EC2	Binocular	SWF10X		Zoom1x-4x		Pillar Big Stand		Frosted Glass			
A23.1201-EC2-3	Trinocular	SWF10X		Zoom1x-4x		Pillar Big Stand		Frosted Glass			
Technical Data											
Eyepiece	Standard Objective 1-4X		Auxiliary Objective								
	0.5X				0.75X		1.5X		2.0X		
	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance	Field of View	Working Distance	
	5X	20-5	91	40-11	170	25-6.5	100	13-3.5	42	10-2.5	28
	10X	23-5.5	91	52-12	170	30-7	100	15-4.8	42	11-5.3	28
	15X	15.5-4	91	36-8.5	170	21-5	100	10.5-2.5	42	8-2	28
20X	10.5-3	91	25-5.8	170	14-3.5	100	7-1.8	42	5.5-1.5	28	
Optional Accessories										Item No.	
Eyepiece	WF5X /22									A51.1201-05	
	WF10X /23									A51.1201-10	
	WF15X/13									A51.1201-15	
	WF20X/10									A51.1201-20	
Auxiliary Objective	0.5X									A52.1282-05	
	0.75X									A52.1282-75	
	1.5X									A52.1282-15	
	2X									A52.1282-20	
Figure Eyepiece	VD-TV									A51.1207-TV	
	VD-PC									A51.1207-PC	
Stage	Darkfield stage,Φ95mm									A54.1202-1	
Plate	White/black Working Plate Φ95mm									A54.1203-2	
Illumination	more choice pls click here:light source										
Packing	Wooden cabinet with lock and key									A50.1202	
	Aluminium case									A50.1203	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China