



## Lab Biological Microscope

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1132
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 400X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope, biological compound microscope**

## Product Description

### Lab Achromatic Biological Microscope Monocular For Student A11.1132

#### Features:

Full metal body, high stability  
Cordless, compact size, easy to carry  
Top and bottom LED light source  
Competitive low price

#### Specification:

A11.1132 Economical Cordless Student Microscope Specification	
Head	Monocular Head
Eyepiece	WF10X
Objective	Achromatic 4X
	Achromatic 10X
	Achromatic 40X, Spring
Nosepiece	Triple Revolving Nosepiece
Stage	Plain Stage with Clips, size 90*90mm
Condenser	6 Hole Disc Diaphragm
Focusing	Coarse Focusing Knobs
Light Source	Top and Bottom LED light, powered by AA 1.5V Battery 3 pcs. (No battery included when ship by air)
Body	Main Body Size 10*15*27cm



A11.1132



Triple Nosepiece, Upper LED Light






6 Hole Disc Diaphragm



Coarse Focusing Knobs



 0086 13911110627    [sale@optoedu.com](mailto:sale@optoedu.com)    [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China