



## Textile Handheld Digital Microscope

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A 34.4185
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Model Number: A34.4185
- Application: Home
- Certification: CE
- Magnification: 200x
- Theory: Handheld Digital Microscope
- Image Sensor: CMOS Sensor
- Highlight: Portable Digital Microscope, hand held microscope



### More Images



Our Product

### Product Description

#### Description:

New Product in 2013!

This model break the limitation of distance, it works within 10 meter's range without any wire!

It is a perfect match with IOS mobile device and android device.

It can be easily used in different fields including:

1. Textile industrial for textile inspection
2. Printing inspection
3. Industrial inspection: PCB, Precision machinery

- 4.Educational Purpose
- 5.Hair Examination
- 6.Skin Examination
- 7.Microbiological Observation
- 8.Jewelry&Coin(Collections)inspection
- 9.Visual Assistance

**Specifications:**

A 34.4185 Specification	
Magnification	5 X ~ 200X ( Manual )
Focus Range	0mm ~40mm,other range can be customized
Image Sensor	CMOS Sensor
Application Software	HotViewer
Resolution	Video Capture Resolution 0.3M Still Image Capture Resolution 640*480,can be magnified and minified
Controller	High Speed DSP(M-jpeg)
Frame Rate	20f/s VGA under 600 lux brightness
Power Supply	Rechargeable 1000mA self-protection lithium battery
Illumination	Built-in 4-8 White-light LED and adjustable illumination ensure the magnified images are clear and bright Scope:0-30000lux
Take Images	Snap shot in build-in application program
Wifi Standard	IEEE802.11 b/g/n Standard compliant
Fixed Foundation	Tailor-made for the microscope
Brightness Control	Manual adjustment
System	Support O/S:WIFI capable such as IOS and Android
Color	Deep Blue
Language	English,other language can be customized
Dimension	142mm(L)*37mm(R)
Net Weight	100g







**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China