



1000x Multi Viewing Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1013-C
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Educational, School
- Certification: CE, Rohs
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Drawtube: Binocular
- Highlight: dual head microscope, Dual View Microscope



More Images



Product Description

Specifications:

A17.1013-C Multi-Viewing Microscope, 5 Position	
Head	5 pcs Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm
Eyepiece	WF10X, 10 pcs
Objective	Achromatic 4x, 10x, 40x(s), 100x(s, oil)
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double Layer Mechanical Stage 140mm x 140mm, Travel Range 75x50mm
Condenser	Abbe N.A.=1.25 with Iris Diaphragm & Filter
Focusing	Coaxial Coarse & Fine Focusing Knobs, Fine division 0.002mm, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.2mm Per Rotation, Moving Range 28mm
Light Source	Base with Built-in Halogen Lamp 24V/100W , with Continue Intensity Control, 110V/220V
Indicator	Demo. Indicator with Illuminator of 6V/3W





A17.1013 Multi-Viewing Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China