



Gems Jewelry Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1203
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Diamond, Gems
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 7.5x - 5x
- Theory: Jewelry Microscope
- Highlight: diamond microscope, gemological microscope



Product Description

Specifications:

A24.1203 Jewelry Microscope			
Head	Inclined Binocular Head,360° Rotatable	●	
	Inclined Trinocular Head,360° Rotatable		●
Eyepiece	SWF10x		
Objective	0.75x~5x		
Incident Light	Halogen Lamp 12V/10W, With Single Fiber		
	Fluorescent 10W Square		
Transmitted Light	Halogen Lamp 12V/10W		
Working Stage	Jewel Tweezer		
	Dark Field Working Stage, With Diaphram		
A24.1203 Jewelry Microscope , Optional Accessories			
Eyepiece	WF5x	E51.1221-05	
	WF10x	E51.1221-10	
	WF15x	E51.1221-15	
	WF20x	E51.1221-20	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China