



## Sliding Binocular Biological Compound Microscope A11.1104

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1104
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Laboratory, Hobby
- Certification: CE
- Drawtube: Binocular, Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,  
binocular biological microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Sliding Binocular Biological Compound Microscope A11.1104

#### Specifications:

A11.1104 Biological microscope		-A	-B	-C	-D
Head	Monocular Head,45° Inclined,360° Rotatable	●	●	●	-
	Sliding Binocular Head,45° Inclined,55mm-75mm	-	-	-	●
Eyepiece	WF10X/18mm	●	●	●	●
Objective	Achromatic 4X	●	●	●	●
	Achromatic 10X	●	●	●	●
	Achromatic 40X, Spring	●	●	●	●
	Achromatic 100X, Spring	-	-	●	●
Nosepiece	Triple Revolving Nosepiece	●	●	-	-
	Quadruple Revolving Nosepiece	-	-	●	●
Stage	Plain Stage with Specimen Clips, Size1210x110mm	●	●	-	-
	Mechanical Stage,Moving Range:60mmx30mm	-	-	●	●
	N.A.=0.65, with Disc Diaphragm	●	-	-	-
Condenser	N.A.=0.65, with Iris Diaphragm	-	●	-	-
	Abbe Condenser N.A.=1.25, with Iris Diaphragm & Filter	-	-	●	●
Focusing	Seperate Coarse & Fine Focusing Knobs	●	●	●	●
Light Source	Incandescent Lamp,220V/20W or 110V/20W	●	●	●	●
<b>Optional Accessories</b>					
Head	Demonstration WP-T				
	Demonstration WP-V				
Eyepiece	WF16X				
	WF20X				
	P16X				
Objective	Achromatic 20X				
	Achromatic 60X(S)				
Light Source	LED				







**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China