



Sliding Binocular Dark Field Microscope Halogen Lamp A11.1110

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1110
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Hobby
- Certification: CE
- Drawtube: Binocular, Monocular, Trinocular
- Illumination: Halogen Lamp
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**



Product Description

Specifications:

A11.1110 Biological Microscope		-A	-B	-C	-D	-E
Head	Monocular Head,45°Inclined, 360°Rotatable	●	-	-	-	-
	Compesation Free Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm	●	-	-	-	-
	Compesation Free Trinocular Head, 30°Inclined, Interpupillary Distance 55-75mm	-	-	●	-	-
	Sliding Binocular Head, 45°Inclined,Interpupillary Distance 55-75mm	-	-	-	●	-
	Sliding Trinocular Head, 45°Inclined,Interpupillary Distance 55-75mm	-	-	-	-	●
Eyepiece	WF10X/18mm	●	●	●	●	●
Objective	Achromatic 4X	●	●	●	●	●
	Achromatic 10X	●	●	●	●	●
	Achromatic 40X, Spring	●	●	●	●	●
	Achromatic 100X, Spring, Oil	●	●	●	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●	●	●	●
Stage	Double Layer Mechanical Stage, Size 140x140mm, Moving Range 70x50mm	●	●	●	●	●
Condenser	Abbe Condenser N.A.=1.25, with Iris Diaphragm and Filter	●	●	●	●	●
Focusing	Coaxial Coarse & Fine Focusing Scale Value 0.002mm	●	●	●	●	●
Illumination	Kohler Illumination	●	●	●	●	●
Light Source	Halogen Lamp 6V 20W,Brightness Adjustable	●	●	●	●	●
Optional Accessories						
Eyepiece	WF16X					
	WF20X					
	P16X					
Objective	Achromatic 20X					
	Achromatic 60X, Spring					
	Semi-plan Objective 4X					
	Semi-plan Objective 10X					
	Semi-plan Objective 20X					
	Semi-plan Objective 40X(S)					
	Semi-plan Objective 60X(S)					
	Semi-plan Objective 100X(S)					
	Plan Objective 4X					
	Plan Objective 10X					
	Plan Objective 20X					
	Plan Objective 40X(S)					
	Plan Objective 60X(S)					
Plan Objective 100X(S)						
Dark Field	Darkfield					
Photo Adapter	Photo Adapter					
CCD Adapter	0.5x					
	0.57x					
	0.75x					



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China