



LED Monocular Biological Microscope With 45°Inclined Head A11.1119

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1119
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Laboratory
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



Product Description

LED Monocular Biological Microscope With 45°Inclined Head A11.1119

Specifications:

A11.1119 Biological Microscope		-A	-B
Head	Monocular Head,45° inclined,360° Rotatable	●	●
Eyepiece	WF10X	●	●
Objective	Achromatic 4X	●	●
	Achromatic 10X	●	●
	Achromatic 40X, Spring	●	●
Nosepiece	Triple Revolving Nosepiece	●	●
Stage	Plain Stage with Specimen Clips, Size 90mmx90mm	●	●
Condenser	Single Lens with Disc Diaphragm	●	●
Focusing	Coarse Focusing	●	-
	Coaxial Coarse & Fine Focusing Adjustment Fine Focusing Scale Value 0.002mm	-	●
Light Source	Incandescent Lamp,220V/20W or 110V/20W	●	●
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
Objective	Achromatic 20X		
	Achromatic 60X(S)		
Light Source	LED		





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China