



## Coaxial Coarse Inverted Optical Microscope Transmitting Light , Rohs A14.2601

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.2601
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Laboratory, Research
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Max Magnification: 400X
- Product Name: Inverted Optical Microscope
- Theory: Inverted Microscope
- Highlight: **inverted biological microscope,**  
**inverted fluorescence microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Specifications:

A14.2602 Inverted Microscope Specification		
Head	Seidentopf Trinocular Head, 45° Inclined, Interpupillary Distance 50-75mm, Light Distribution (binocular tube:photo tube): Both 100:0 and 20:80	
Eyepiece	WF10x/20mm	
	Telescope	
Long Working Distance Objective	Plan Objective LPL10x/0.25 WD=7.9mm	
	Plan Objective LPL25x/0.40 WD=5.0mm	
	Plan Objective LPL40x/0.65 WD=3.0mm	
	Plan Phase Contrast Objective LPL Ph+ 25x/0.40 WD=5.0mm	
Focus System	Low positioned, coaxial coarse & fine focusing system, coarse 42mm/rot, fine 0.2mm/rot, minimum reading 2µm	
Nosepiece	Quadruple	
Stage	Mechanical Stage 206(X) x 170(Y)mm, with attached mechanical stage moving range 128(X) x 80(Y)mm	
Culture Dish Holder	65mm 35mm	
Condenser	Slide Relief Phase Contrast Condenser N.A. 0.35	
Phase Contrast	Phase Contrast Slider 25x	
Filter	Blue, Green & Gray	
Light Source	Transmitting light, Kohler Illumination System: 12V/20W Halogen	
Power Supply	85-265V, 50/60Hz	
Optional Accessories	Item No.	
Eyepiece	WF16x/13mm	A51.2602-1613A
	10x Reticle Eyepiece (0.1mm)	A51.2605-10
Long Working Distance Objective	Plan Objective LPL2.5x/0.07, WD=29.7mm	A52.2608-2.5
	Plan Objective LPL4x/0.10, WD=17.5mm	A52.2608-4
	Plan Phase Contrast Objective LPL Ph+ 10x/0.25 W.D.=7.9 mm	A5C.2630-10
	Plan Phase Contrast Objective LPL Ph+ 40X/0.25 W.D.=3.0 mm	A5C.2630-40
Phase Contrast	Phase Contrast Slider 10x	A5C.2611-10
	Phase Contrast Slider 40x	A5C.2611-40
Culture Dish Holder	54mm	A54.2606-54
Stage	Attached Mechanical Stage for 96 holes Plate	A54.2601
Adapter	C-Mount Adapter 0.6x (Focus Adjustable)	A55.2601-B06
	C-Mount Adapter 1x	A55.2601-B1

Digital SLR Camera Adapter	For Nikon	A55.2622-N
	For Sony	A55.2622-S
	For Canon	A55.2622-C
Light Source	Upgrade to 6V30W Halogen	A56.2650-30W
Other	Micrometer 0.1mm	A58.2601-M






**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China