



40x - 600x Metal Polarizing Light Microscope With Rotary Stage A15.1103

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1103
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Metal
- Certification: CE, Rohs
- Drawtube: Binocular
- Magnification: 40x - 600x
- Theory: Polarizing Light Microscope
- Stage: Rotary Stage
- Highlight: **student compound microscope ,
highest magnification microscope**



Product Description

40x - 600x Metal Polarizing Light Microscope With Rotary Stage A15.1103

Specifications:

A15.1103 Polarizing Microscope Specification		-A	-B
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●
Eyepiece	Wide field WF10X/22mm	●	●
	Gridding Eyepiece	●	●
	Cross Eyepiece	●	●
	Eyepiece with Reticule	●	●
Objective	Stress-Free Infinity 40X	●	●
	Plan-Achromatic		
	Metallurgical	●	●
	Objective		
	4X	●	●
	Stress-Free Infinity 10X	●	●
	Plan Achromatic	●	●
	Objective	●	●
Condenser	Swing Type N.A.1.25 Abbe Condenser with Iris Diaphragm	●	●
	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm.	●	●
Nosepiece	Quintuple Nosepiece with centering adjustable	●	●
Filter	Green Inserting Type Filting	●	●
	Blue Inserting Type Filting	●	●
	Neutral Inserting Type Filting	●	●
Stage	Rotary Stage,Stage Size: Φ172mm,Moving Range:360°	●	●
Testing Piece	Gypsum 1 λTesting Piece	●	●
	Mica1/4 λ	●	●
	Quartz Wedge Testing Piece	●	●
Illumination	High brightness Fixed Illumination	●	●
Light Source	Upper:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness	●	●
	Down:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness	●	●
Checking Tool	0.01mm Micrometer	●	●
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
	P16X		
Stress-Free Infinity Plan-Achromatic Metallurgical Objective	4X		
	10X		
	20X		
Photo Adapter	Photo Adapter		
CCD Adapter	0.5X		
	0.57X		
	0.75X		



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China