



## 40x - 600x Metal Polarizing Light Microscope With Rotary Stage A15.1103

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1103
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Application: Metal
- Certification: CE, Rohs
- Drawtube: Binocular
- Magnification: 40x - 600x
- Theory: Polarizing Light Microscope
- Stage: Rotary Stage
- Highlight: **student compound microscope ,  
highest magnification microscope**



## Product Description

### 40x - 600x Metal Polarizing Light Microscope With Rotary Stage A15.1103

#### Specifications:

A15.1103 Polarizing Microscope Specification		-A	-B		
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-		
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●		
Eyepiece	Wide field WF10X/22mm	●	●		
	Gridding Eyepiece	●	●		
	Cross Eyepiece	●	●		
	Eyepiece with Reticule	●	●		
Objective	Stress-Free Infinity	40X	●	●	
		Plan-Achromatic Metallurgical	60X	●	●
	Objective	4X	●	●	
		Stress-Free Infinity	10X	●	●
		Plan Achromatic	20X	●	●
		Objective	40X	●	●
		Objective	60X	●	●
		Condenser	Swing Type N.A.1.25 Abbe Condenser with Iris Diaphragm	●	●
Focus System	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm.	●	●		
Nosepiece	Quintuple Nosepiece with centering adjustable	●	●		
Filter	Green Inserting Type Filting	●	●		
	Blue Inserting Type Filting	●	●		
	Neutral Inserting Type Filting	●	●		
Stage	Rotary Stage,Stage Size: Φ172mm,Moving Range:360°	●	●		
Testing Piece	Gypsum 1 λ Testing Piece	●	●		
	Mica1/4 λ Quartz Wedge Testing Piece	●	●		
Illumination	High brightness Fixed Illumination	●	●		
Light Source	Upper:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness	●	●		
	Down:Halogen Bulb 12V/30W,AC 85V-230V ,Adjustble Brightness	●	●		
Checking Tool	0.01mm Micrometer	●	●		
Optional Accessories					
Eyepiece	WF16X				
	WF20X				
	P16X				
Stress-Free Infinity Plan-Achromatic Metallurgical Objective	4X				
	10X				
	20X				
Photo Adapter	Photo Adapter				
CCD Adapter	0.5X				
	0.57X				
	0.75X				



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China