



Research Inverted Phase Contrast Microscope 400X With Coaxial Coarse A19.2602

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A19.2602
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular
- Theory: Phase Contrast Microscope
- Nosepiece: Quadruple
- Light Source: Transmitting Light, Kohler Illumination System: 12V/20W Halogen
- Highlight: **phase contrast light microscope, phase contrast microscopes**



Product Description

Research Inverted Phase Contrast Microscope With Coaxial Coarse A19.2602

Specifications:

A19.2602 Inverted Phase Contrast Microscope Specification		
Head	Seidentopf Trinocular Head, 45° Inclined, Interpupillary Distance 50-75mm, Light Distribution (binocular tube: photo tube): Both 100:0 and 20:80	
Eyepiece	WF10x/20mm	
	WF16x/16mm	
	10x Reticle Eyepiece (0.1mm)	
	Telescope	
Long Working Distance Objective	Plan Objective LPL4x/0.10 WD=17.5mm	
	Plan Objective LPL25x/0.40 WD=5.0mm	
	Plan Objective LPL40x/0.65 WD=3.0mm	
	Plan Phase Contrast Objective LPL Ph+ 10x/0.25 WD=7.9mm	
	Plan Phase Contrast Objective LPL Ph+ 25x/0.40 WD=5.0mm	
	Plan Phase Contrast Objective LPL Ph+ 40x/0.65 WD=3.0mm	
Focus System	Low positioned, coaxial coarse & fine focusing system, coarse 42mm/rot, fine 0.2mm/rot, minimum reading 2μm	
Nosepiece	Quadruple	
Stage	Mechanical Stage 206(X) x 170(Y)mm, with attached mechanical stage moving range 128(X) x 80(Y)mm	
Culture Dish Holder	65mm	
	35mm	
Condenser	Turnable Relief Phase Contrast Condenser N.A.0.35 W.D. =45mm	
Phase Contrast	Phase Contrast Slider 10x	
	Phase Contrast Slider 25x	
	Phase Contrast Slider 40x	
Filter	Blue, Green & Gray	
Light Source	Transmitting light, Kohler Illumination System: 12V/20W Halogen	
Adapter	C-Mount Adapter 1x	
Power Supply	85-265V, 50/60Hz	
Optional Accessories		Item No.
Long Working Distance Objective	Plan Objective LPL2.5x/0.07, WD=29.7mm	A52.2608-2.5
	Plan Objective LPL4x/0.10, WD=17.5mm	A52.2608-4
	Plan Phase Contrast Objective LPL Ph+ 10x/0.25 W.D.=7.9 mm	A5C.2630-10
	Plan Phase Contrast Objective LPL Ph+ 40X/0.25 W.D.=3.0 mm	A5C.2630-40
Culture Dish Holder	54mm	A54.2606-54
Stage	Attached Mechanical Stage for 96 holes Plate	A54.2601
Adapter	C-Mount Adapter 0.6x (Focus Adjustable)	A55.2601-B06
Digital SLR Camera Adapter	For Nikon	A55.2622-N
	For Sony	A55.2622-S
	For Canon	A55.2622-C
Light Source	Upgrade to 6V30W Halogen	A56.2650-30W
Other	Micrometer 0.1mm	A58.2601-M



F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China