



## Binocular / Trinocular Stereo Optical Microscope For Jewelry A22.1001

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1001
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Gem, Jewelry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,  
stereo dissecting microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Binocular / Trinocular Stereo Optical Microscope For Jewelry A22.1001

#### Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

#### Specifications:

A22.1001 Stereo Microscope Microscope		13	14	12	23	24
Head	Binocular head,inclinedat 45°	•	•	•	•	•
Eyepiece	EW10X /23mm	•	•	•	•	•
Objective	1x/3x	•				
	1x/4x		•			
	1x/2x			•		
	2x/3x				•	
	2x/4x					•
Working Dist ance	100mm	•	•	•	•	•
Illumination	Stand with LEDtransmission & reflectiong illumination	•	•	•	•	•
A22.1001 Stereo Microscope Microscope Optional Accessories						
Head	Binocular head,inclined at 60°	A22.1001-BH				
	Trinocular head,inclined at 45°	A22.1001-TH				
Eyepiece	WF15X/16mm	A51.1029-15				
	WF30X	A51.1029-30				
	High-point WF20X/12mm	A51.1029-20				
	Low-point WF20X/12mm	A51.1029-21				
	WF10XEyepiece with cross hair	A51.1029-10				
Auxiliary objective	2X	A52.1026-2				
	0.75X	A52.1026-75				
	0.63X	A52.1026-63				
	0.5X	A52.1026-5				
Stand	Stand without illumination	A54.1022-1				
	Stand with LED transmission illumination	A54.1022-2				
	Stand with LED reflectiong illumination	A54.1022-3				
Video Attach ment	C-Mount 0.5x	A55.1002-05				
	C-Mount 1.0x	A55.1002-10				

#### About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China