



Stereo Optical Microscope Stereo Zoom Microscopes Halogen Lamp A22.1102

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Medical|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp|LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo dissecting microscope**



Product Description

Halogen Lamp Stereo Optical Microscope Stereo Zoom Microscopes A22.1102

Description:

- WF10X/20 eyepiece (standard)
- Top and Bottom light 12V/15W halogen Lamp, brightness is adjustable
- LED illumination or Rechargeable LED illumination for option
- WF15X/15 WF20X/11 for option

Specifications:

Model No.	Objective	Eyepiece				Working Distance (mm)	
		WF10X/20(standard)	WF15X/15(Optional)	WF20X/11(Optional)			
		Mag.	FOV (mm)	Mag.	FOV (mm)	Mag.	FOV (mm)
A22.1102 2X		20X	11.5	30X	7.5	40X	6.75
-C	4X	40X	5.75	60X	3.75	80X	3.37
A22.1102							110
-A							
A22.1102							
-B							



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China