



Training Zoom Stereo Optical Microscope 7x - 45x , Rohs A23.0901-S1

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0901-S1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x - 45x
- Applications: Industry, Jewelry
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: stereo zoom microscopes,
stereo dissecting microscope



Product Description

Training Zoom Stereo Optical Microscope 7x - 45x , Rohs A23.0901-S1

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Single Arm Universal Pole Stand

Specifications:

A23.0901-S1 Zoom Stereo Microscope Specification							B	T
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
Eyepiece	WF10x/20mm, High Eyepoint							
Zoom Lens	0.7~4.5x							
Zoom Ratio	1:6.5							
Magnification	7x-45x, Up to 3.5x-180x With Optional Eyepieces and Auxiliary Objectives							
Working Distance	100mm, With Auxiliary Objectives 30~165mm							
Field of View	1~95.2mm							
Light Source	No Light Source							
Stand	Single Arm Universal Pole Stand STL1, Height 308mm, Focusing Pole Height 104mm, Level Move Range 235mm, Base 230*230*40mm,							
A23.0901-S1 Zoom Stereo Microscope Optical Data Sheet								
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W.D.	100mm		165mm		45mm		30mm	
Eyepiece	Mag	View field	Mag	View field	Mag	View field	Mag	View field
	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm
10X/20mm	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm
	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm
15X/15mm	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm
	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm
20X/10mm	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm
A23.0901-S1 Zoom Stereo Microscope Optional Accessories							Cata. No.	
Head	Binocular Head + WF10x/20mm						A23.0901-BH	
	Trinocular Head + WF10x/20mm						A23.0901-TH	
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm						A51.0922-10	
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm						A51.0922-15	
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm						A51.0922-20	
	WF10x/20mm,with Reticle						A51.0932-10	
	WF15x/15mm,with Reticle						A51.0932-15	
	WF20x/10mm,with Reticle						A51.0932-20	
Auxiliary Lens	0.5X,W.D. 165mm						A52.0927-05	
	0.75X,W.D. 105mm						A52.0927-75	
	1.5X,W.D. 45mm						A52.0927-15	
	2.0X,W.D. 30mm						A52.0927-20	
CCD Adapter	0.50x, For 1/2" CCD						A55.0932-05	
	0.35x, For 1/3" CCD						A55.0932-35	



A23.0901-S1 Zoom Stereo Microscope 0.7~4.5x




Eyepiece WF10x/20mm, High Eyepoint



Magnification Zoom Knob 0.7~4.5x



Trinocular Head



Universal Stand STL1



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China