



Industry Learning Stereo Zoom Microscopes For Coin / Stamp A23.0901-S2

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0901-S2
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x - 45x
- Applications: Industry, Jewelry
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes, stereo dissecting microscope**



Product Description

Industry Learning Stereo Zoom Microscopes For Coin / Stamp A23.0901-S2

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Double Arm Universal Pole Stand

Specifications:

A23.0901-S2 Zoom Stereo Microscope Specification							B	T
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
Eyepiece	WF10x/20mm, High Eyepoint							
Zoom Lens	0.7~4.5x							
Zoom Ratio	1:6.5							
Magnification	7x-45x, Up to 3.5x-180x With Optional Eyepieces and Auxiliary Objectives							
Working Distance	100mm, With Auxiliary Objectives 30~165mm							
Field of View	1~95.2mm							
Light Source	No Light Source							
Stand	Double Arm Universal Pole Stand STL2, Height 380mm, Focusing Pole Height 140mm, Level Move Range 235mm, Base 210*260*20mm,							
A23.0901-S2 Zoom Stereo Microscope Optical Data Sheet								
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X	
W.D.	100mm		165mm		45mm		30mm	
Eyepiece	Mag	View field	Mag	View field	Mag	View field	Mag	View field
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm
A23.0901-S2 Zoom Stereo Microscope Optional Accessories							Cata. No.	
Head	Binocular Head + WF10x/20mm						A23.0901-BH	
	Trinocular Head + WF10x/20mm						A23.0901-TH	
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm						A51.0922-10	
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm						A51.0922-15	
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm						A51.0922-20	
	WF10x/20mm,with Reticle						A51.0932-10	
	WF15x/15mm,with Reticle						A51.0932-15	
Auxiliary Lens	WF20x/10mm,with Reticle						A51.0932-20	
	0.5X,W.D. 165mm						A52.0927-05	
	0.75X,W.D. 105mm						A52.0927-75	
	1.5X,W.D. 45mm						A52.0927-15	
CCD Adapter	2.0X,W.D. 30mm						A52.0927-20	
	0.50x, For 1/2" CCD						A55.0932-05	
	0.35x, For 1/3" CCD						A55.0932-35	



Trinocular Head



Binocular Head
WF10x/20mm Eyepiece, High Eyepoint



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China