



Medical Stereo Optical Microscope Fluorescent Light Microscopes CE A23.0901-ST

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0901-ST
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x - 45x
- Applications: Jewelry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**



More Images



Product Description

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Up 12V 15W Halogen & Bottom 230V 7W Fluorescent Light
- 110~240V Wide Range Voltage

Specifications:

A23.0901-ST Zoom Stereo Microscope Specification						2B	2T	3B	3T
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm					o	o		
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm					o	o		
Eyepiece	WF10x/20mm, High Eyepoint								
Zoom Lens	0.7~4.5x								
Zoom Ratio	1:6.5								
Magnification	7x-45x, Up to 3.5x-180x With Optional Eyepieces and Auxiliary Objectives								
Working Distance	100mm, With Auxiliary Objectives 30~165mm								
Field of View	1~95.2mm								
Stand	Pole Stand ST2/ST3, Height 248mm, Base 200*255*60mm								
Light Source	Up 12V 15W Halogen & Bottom 12V 15W Halgen Light					o			
	Up 12V 15W Halogen & Bottom 230V 7W Fluorescent Light					o			
Power	110~240V Wide Range Voltage								
A23.0901-ST Zoom Stereo Microscope Optical Data Sheet									
Objective	Standard 1.0x		Auxiliary 0.5X		Auxiliary 1.5X		Auxiliary 2X		
W.D.	100mm		165mm		45mm		30mm		
Eyepiece	Mag	View field	Mag	View field	Mag	View field	Mag	View field	
10X/20mm	7.0X	28.6mm	3.5X	57.1mm	10.5X	19mm	14X	14.3mm	
	45.0X	4.4mm	22.5X	8.9mm	67.5X	3mm	90X	2.2mm	
15X/15mm	10.5X	21.1mm	5.25X	42.8mm	15.75X	14.3mm	21X	10.7mm	
	67.5X	3.3mm	33.75X	6.7mm	101.25X	2.2mm	135X	1.7mm	
20X/10mm	14.0X	14.3mm	7.0X	28.6mm	21.0X	9.5mm	28X	7.1mm	
	90.0X	2.2mm	45.0X	4.4mm	135.0X	1.5mm	180X	1.1mm	
A23.0901-ST Zoom Stereo Microscope Optional Accessories									
							Cata. No.		
Head	Binocular Head + WF10x/20mm						A23.0901-BH		
	Trinocular Head + WF10x/20mm						A23.0901-TH		
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm						A51.0922-10		
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm						A51.0922-15		
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm						A51.0922-20		
	WF10x/20mm,with Reticle						A51.0932-10		
	WF15x/15mm,with Reticle						A51.0932-15		
	WF20x/10mm,with Reticle						A51.0932-20		
Auxiliary Lens	0.5X,W.D. 165mm						A52.0927-05		
	0.75X,W.D. 105mm						A52.0927-75		
	1.5X,W.D. 45mm						A52.0927-15		
	2.0X,W.D. 30mm						A52.0927-20		
CCD Adapter	0.50x, For 1/2" CCD						A55.0932-05		
	0.35x, For 1/3" CCD						A55.0932-35		



Trinocular Head



Binocular Head
WF10x/20mm Eyepiece, High Eyepoint



A23.0901-ST Zoom Stereo Microscope
Up & Bottom Light Source, Brightness Continuously Adjustable



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China