



Clinic Stereo Optical Microscope Stereo Zoom Microscope A23.0902-S1

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0902-S1
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x - 45x
- Applications: Jewelry|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular|Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



Product Description

Clinic Stereo Optical Microscope Stereo Zoom Microscope A23.0902-S1

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Single Arm Universal Pole Stand

Specifications:

A23.0902-S1 Zoom Stereo Microscope Specification							B1	T1
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm						o	
Eyepiece		WF10x/22mm, High Eyepoint						
Zoom Lens		0.67~4.5x						
Zoom Ratio		1:6.7						
Magnification		6.7x~45x, (2x~225x with Auxiliary Objective)						
Working Distance		100mm, With Auxiliary Objectives 26~287mm						
Light Source		No Light Source						
Stand		Single Arm Universal Pole Stand STL1, Height 308mm, Focusing Pole Height 104mm, Level Move Range 235mm, Base 230*230*40mm,						
A23.0902-S1 Zoom Stereo Microscope Optical Data Sheet								
Auxiliary Lens		0.3x	0.4x	0.5x	0.75x	1.0x	1.5x	2.0x
Working Distance		287mm	217mm	177mm	117mm	100mm	47mm	26mm
WF 10x	Magnification	2.0-13.5x	2.7-18.0x	3.4-22.5x	5.0-33.8x	6.7-45.0x	10.1-67.5x	13.4-90.0x
	View Field(mm)	109.5-16.3	82.1-12.2	65.7-9.8	43.8-6.5	32.8-4.9	21.9-3.3	16.4-2.4
WF 15x	Magnification	3.0-20.3x	4.0-27.0x	5.0-33.8x	7.5-50.6x	10.1-67.5x	15.1-101.3x	20.1-135x
	View Field(mm)	79.6-11.9	59.7-8.9	47.8-7.1	31.8-4.7	23.9-3.6	15.9-2.4	11.9-1.8
WF 20x	Magnification	4.0-27.0x	5.4-36.0x	6.7-45.0x	10.1-67.5x	13.4-90.0x	20.1-135.0x	26.8-180x
	View Field(mm)	59.7-8.9	44.8-6.7	35.8-5.3	12.9-2.4	17.9-2.7	11.9-1.8	9-1.3
WF 25x	Magnification	5.0-33.8x	6.7-45.0x	8.4-56.3x	12.6-84.4x	16.8-112.5x	25.1-168.8x	33.5-225x
	View Field(mm)	44.8-6.7	33.6-5.0	26.9-4.0	17.9-2.7	13.4-2.0	9.0-1.3	6.7-1.0
A23.0902-S1 Zoom Stereo Microscope Optional Accessories							Cata. No.	
Head	Binocular Head + WF10x/20mm						A23.0902-BH	
	Trinocular Head + WF10x/20mm						A23.0902-TH	
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm						A51.0923-10	
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm						A51.0923-15	
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm						A51.0923-20	
	WF10x/20mm,with Reticle						A51.0933-10	
	WF15x/15mm,with Reticle						A51.0933-15	
Auxiliary Lens	WF20x/10mm,with Reticle						A51.0933-20	
	0.5X,W.D. 165mm						A52.0928-05	
	0.75X,W.D. 105mm						A52.0928-75	
	1.5X,W.D. 45mm						A52.0928-15	
CCD Adapter	2.0X,W.D. 30mm						A52.0928-20	
	0.50x, For 1/2" CCD						A55.0933-1/2	
0.35x, For 1/3" CCD						A55.0933-1/3		



Trinocular Head



Binocular Head
WF10x/22mm Eyepiece, High Eyepoint



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China